



**Das Deutsche Zentrum für
Neurodegenerative Erkrankungen**

**Introduction to Bioimage analysis:
Retreat 2022 Tuebingen (Constance)**

31/05/2022


Miguel Fernandes

**DZNE e. V. – Bonn / CRFS - Image and Data Analysis
Facility (IDAF)**

Survey


What is your level of experience with analysis?

11 respostas




From this list which software/tools have you used? Select all that apply

11 respostas




Which kind of data are you collecting? Select all that apply

11 respostas



Which modalities of imaging are you using?

11 respostas



Talk outline:

- What is bioimage analysis?
- Basics of image processing
- FIJI: macros/ Napari
- What can we improve?
 - GPU-accelerated image analysis (CLIJ2)
 - AI vs ML vs DL
 - Trainable Weka Segmentation
 - CellProfiler and Cell Analyst
 - Visualization of results and statistical analysis with Python
- Do it yourself with your data!

What is bioimage analysis?

One possible definition:

The extraction of relevant information from digital images, at large scale and high throughput in the context of biological research

- Precise
- Unbiased
- Reproducible
- Scalable


```
31
32 self.file = None
33 self.fingerprints = set()
34 self.logdups = True
35 self.debug = debug
36 self.logger = logging.getLogger(__name__)
37 if path:
38     self.file = open(os.path.join(path, "requests.log"),
39                     "a")
40     self.fingerprints.update([request for r in self.requests])
41
42 @classmethod
43 def from_settings(cls, settings):
44     debug = settings.getbool("debug")
45     return cls(job_dir(settings), debug)
46
47 def request_seen(self, request):
48     fp = self.request_fingerprint(request)
49     if fp in self.fingerprints:
50         return True
51     self.fingerprints.add(fp)
52     if self.file:
53         self.file.write(fp + os.linesep)
54
55 def request_fingerprint(self, request):
56     return request_fingerprint(request)
```

Image from Unsplash

What is the goal of Bioimage analysis?

- *Obtain useful insights out of complicated and heterogeneous image and related metadata.*

See https://en.wikipedia.org/wiki/Bioimage_informatics for details)



From <https://www.nature.com/articles/nbt.3722>

Meijering, E. et al. "Imaging the future of bioimage analysis."

Computational analysis is becoming the bottleneck

Seeing Is Believing: Quantifying Is Convincing: Computational Image Analysis in Biology


Ivo F Sbalzarini ^{1 2 3}

Affiliations + expand

PMID: 27207361 DOI: [10.1007/978-3-319-28549-8_1](https://doi.org/10.1007/978-3-319-28549-8_1)


Abstract

Imaging is center stage in biology. Advances in microscopy and labeling techniques have enabled unprecedented observations and continue to inspire new developments. Efficient and accurate quantification and computational analysis of the acquired images, however, are becoming the bottleneck. We review different paradigms of computational image analysis for intracellular, single-cell, and tissue-level imaging, providing pointers to the specialized literature and listing available software tools. We place particular emphasis on clear categorization of image-analysis frameworks and on identifying current trends and challenges in the field. We further outline some of the methodological advances that are required in order to use images as quantitative scientific measurements.



<https://pubmed.ncbi.nlm.nih.gov/27207361/>

Example of data analysis



Example Cellpose workflow

Many different tools and steps

How to approach?

Many resources freely available:

<https://bioimagebook.github.io/>


<https://forum.image.sc/>

<http://eubias.org/NEUBIAS/>

<http://biii.eu/>



<https://www.youtube.com/c/haesleinhuepf/videos>

... and more



Modified from https://github.com/marionlouveau/NEUBIAS_TS15_DataViz


Example today: Image J / FIJI



Fiji (Fiji Is Just ImageJ):
distribution of ImageJ2, bundling a lot of
plugins which facilitate scientific image
analysis

Example CLIJ2 analysis workflow
<https://github.com/clij/clij2>


Alternative to FIJI: Napari (Python)



napari is a fast, interactive, multi-dimensional image viewer for Python

<https://napari.org/>

Napari (Python)



You can write scripts as well for automation

```
from skimage.io import imread
import pyclesperanto_prototype as cle
import napari_simpleitk_image_processing as nsitk
import napari

if 'viewer' not in globals():
    viewer = napari.Viewer()

image0_m = viewer.layers['membrane'].data




# gaussian blur
image1_gb = cle.gaussian_blur(image0_m, None, 1.0, 1.0, 0.0)

viewer.add_image(image1_gb, name='Result of gaussian_blur
(clesperanto)')
```

<https://github.com/BiAPoL/napari-clusters-plotter>

<https://napari.org/>

Installing FIJI





You can even run ImageJ on your browser

<https://ij.imjoy.io/>

From <https://imagej.net>

Updating FIJI and installing plugins, scripts and macros





From <https://imagej.net>

Check properties of the image using FIJI

Image > Properties

and

Image > Show Information



This is a composite confocal image of primary hippocampal neurons. High affinity bungarotoxin receptors were stained with fluorescent bungarotoxin (c=1). The nicotinic acetylcholine alpha7 subunit is immunofluorescently labeled (c=2). A nAChR chaperone protein fused with CFP was transiently transfected into the neurons (c=3). Nuclei were dyed with Hoechst (c=4). A Nomarski optics image shows the morphology of the neuron (c=5). Image is courtesy of John Alexander.

(Fiji is Just) ImageJ 2.1.0/1.53c; Java 1.8.0_202 [64-bit]; Mac OS X 10.16; 877MB of 14766MB (5%)

Title: Rat_Hippocampal_Neuron.tif
Width: 81.92 µm (512)
Height: 81.92 µm (512)
Size: 2.5MB
Resolution: 6.25 pixels per µm
Voxel size: 0.16x0.16x1 µm³
ID: -5
Bits per pixel: 16 (unsigned)
Display ranges
1: 472-2436
2: 548.3125-2935.75
3: 504.8477-942.6484
4: 518.3594-3141.3472
5: 1937.9375-3136.4940
Image: 1/5
Channels: 5
Composite mode: "composite"
No threshold
ScaleToFit: false
Uncalibrated
URL: http://imagej.nih.gov/ij/images/Rat_Hippocampal_Neuron.zip
Screen location: 429,389 (3440x1440)
SetMenuBarCount: 18 (0ms)
Coordinate origin: 0,0,0
No properties
No overlay
No selection

Basics of image processing


What is a pixel?

A pixel (“picture element”) is the smallest unit of a digital image or graphic that can be displayed and represented on a digital display device.

Pixels are combined to form a complete image, video, text, or any visible thing on a computer display.

A Pixel Is Not A Little Square!!!

A pixel is a data point, a point sample. A pixel does not need to be rendered as a small square.



<https://en.wikipedia.org/wiki/Pixel>


<https://bioimagebook.github.io/>

Pixel size vs resolution

Pixel size is a property you configure during imaging at the microscope

Resolution is a property of your imaging system.

Spatial resolution: measure of how closely lines can be located in an image while still being differentiable




<https://bioimagebook.github.io>

Voxels (3D data)

Voxel: “Volume element”, the smallest unit of the sampled 3D volume

Usually anisotropic: not the same in all directions

The section thickness together with the xy-pixel dimension defines the voxel size.



Bit depth

The bit-depth n enumerates how many different intensity values are present in an image: 2^n grey values

In microscopy, images are usually stored as 8, 12 or 16-bit images.

Higher bit depths corresponding to a greater range of useful image information available





Image from <http://zeiss-campus.magnet.fsu.edu/articles/basics/digitalimaging.html>

Examples of preprocessing the images before further analysis



Example mean filter



Implementing a 3x3 mean filter



Original image




<https://bioimagebook.github.io/>

Removing noise / preprocessing the data before segmentation is important



Today we are going to have a hands-on on segmentation and labelling

Segmentation and labelling





Macro recorder

Create our first macro

The recorder

Typically, macros are not written from scratch, but recorded using the Macro Recorder: Just click on *Plugins > Macros > Record...* and perform some actions. These actions will be recorded in the recorder window, and you can click **Create** to open the recorded instructions in an editor:



Example script is provided
(*macro_connected_component_v1.ijm*)

From <https://imagej.net>

What if you want to apply these operations to all images in a given directory?

```
1 /* This macro counts the number of objects
2 * in a folder with multiple 2D images
3 *
4 Author: Miguel Fernandes (IDAF)
5 Date: 24 May 2022
6 */
7 /*
8 Images
9 There should not be any " " in your directory or file names
10 At the moment working for 2D image data
11
12 Only one channel is expected to perform analysis
13 Keep it consistent across experiments!
14
15 Run this macro and select a folder containing imaging files
16 (inputFolder)
17
18 Select type of input stacks with the format defined in suffix
19 (example .tif files)
20
21 Select output folder to save results for each individual file and summary across all images
22 (outputFolder)
23
24 At the moment only for 2D but in principle could be extended to 3D
25
26 TODO: Improve documentation for each step
27
28 */
29
30 File label="Select a folder to process", style="Directory" inputFolder
31 String label = "File suffix input folder", value = ".tif" suffix
32 File label="Select a folder to save results", style="Directory" outputFolder
33
34
35
36 runFile.getDirectory(inputFolder);
37
38
39 print("User selected input folder: " + inputFolder);
40 print("User selected save folder: " + outputFolder);
41
42
43 setBatchMode(true); //Batch mode on
44
```

Example script is provided
(*macro_connected_component_
v1_BATCH.ijm*)

Now is time to perform the analysis

FIJI time!!!



Your samples?

15 minutes. Lets try it out!


What needs to be improved? Ideas?

- Improve preprocessing and thresholds: what would you try?
- Accelerate imaging processing (e.g. GPU)
- Use Machine learning (ML) to classify our ROIs (Weka)
- Automated analysis for high-content assays (combine classical AI + ML). CellProfiler/Cell Analyst
- Plotting the results and perform statistical tests (Python, see Jupyter Notebook)

GPU Accelerates your image analysis

For installation CLIJ2

<https://clij.github.io/clij2-docs/installationInFiji>



```

1 // To make this script run in Fiji, please activate
2 // the clij and clij2 update sites in your Fiji
3 // installation. Read more: https://clij.github.io
4
5 //Following an update site:
6 //https://imagej.net/update-sites/following#Add update sites
7
8
9 // Generator version: 2.5.1.1
10
11 // Init GPU
12 run("CLIJ2 Macro Extensions", "cl_device=");
13
14 // Load sample dataset Neuron
15 run("Neuron (5 channels)");
16
17 image_1 = getTitle();
18 Ext.CLIJ2_pushCurrentZStack(image_1);
19 // The following auto-generated workflow is made for processing a 2D or 3D dataset.
20 // For processing multiple channels or time points, you need to program a for-loop.
21 // You can learn how to do this online: https://www.youtube.com/watch?v=ulSq-x5\_in4
22
23 // Copy
24 Ext.CLIJ2_copy(image_1, image_2);
25 Ext.CLIJ2_release(image_1);
26
27 Ext.CLIJ2_pull(image_2);
28


```

Example script is provided
(*macro_Neuron_demo_clij2.ijm*)

Haase et al 2020

Licence CC BY 4.0: <https://creativecommons.org/licenses/by/4.0/legalcode>


Artificial intelligence vs Machine learning vs Deep learning




Classical AI vs ML

At the high-level:

- a) either you create rules (Classical AI)
- or
- b) you give examples with labels (Machine learning)



Trainable Weka Segmentation





Lets try it out!

15 minutes

<https://imagej.net/plugins/tws/>

CellProfiler and CellProfiler Analyst

Lets try it out!
15 minutes




<https://cellprofiler.org/>

CP4_MF_workshop.cpproj and
DefaultDB_MF_workshop.workspace provided

<https://cellprofileranalyst.org/>

Visualization of results and statistical analysis (Python and Jupyter Notebook)

Pandas, Matplotlib, Seaborn



Pingouin is an open-source statistical package written in Python 3

<https://pingouin-stats.org/#>

```
pg.ttest(df_1h.Area, df_4h.Area)
```

| | T | dof | tail | p-val | CI95% | cohen-d | BF10 | power |
|--------|----------|------------|-----------|--------------|-------------|----------|-----------|----------|
| T-test | 4.995508 | 669.592392 | two-sided | 7.496529e-07 | [0.0, 0.01] | 0.341099 | 1.389e+04 | 0.998789 |

Is the data normal?

```
print(pg.normality(df_1h.Area))  
print(pg.normality(df_4h.Area))
```

| | | | |
|------|----------|--------------|--------|
| | W | pval | normal |
| Area | 0.976816 | 0.000002 | False |
| | W | pval | normal |
| Area | 0.95573 | 4.555235e-10 | False |

Pingouin can be
used for statistical
analysis

10 minutes

Notebook with analysis is provided

Thank you!!!

Lets use the remaining time to
perform analysis on your own
datasets/discuss about what you
could use for your project

Additional Resources

- https://git.mpi-cbg.de/rhaase/lecture_applied_bioimage_analysis_2020
- <https://montpellierresourcesimagerie.github.io/mri-workshop-machine-learning/>
- <https://imagej.nih.gov/ij/docs/examples/index.html>
- <https://petebankhead.gitbooks.io/imagej-intro/content/>
- Bioimage Data Analysis Workflows book: <https://www.springer.com/gp/book/9783030223854>
- See also Bioimage Analysis: Recommended Reading and Viewing: https://github.com/amgfernandes/Workshop_May_2022_Tuebingen/blob/main/README.md